



Be your best. Get some rest.

☎ 1300 732 695
✉ bocsleepcare@boc.com
☎ 1300 303 253
💻 www.bocsleepcare.com.au

© The Linde Group 2019. Reproduction without permission is strictly prohibited.

Details given in this document are believed to be correct at the time of printing. Whilst proper care has been taken in the preparation, no liability for injury or damage resulting from its improper use can be accepted. Always read and refer to the user manual before operating equipment.



Printed on carbon neutral recycled paper.

MP18-0123 EQA 0619 V1

Wake up to a new way of sleeping

Quality sleep therapy, made affordable



→ Be your best. Get some rest.

Introducing Linde Noctivance®

- ✓ Quiet
- ✓ Efficient
- ✓ Affordable

Perfecting your sleep, breath by breath

Linde Noctivance® APAP (Automatic Positive Airway Pressure) device quietly adjusts pressure to your breathing needs to prevent your airways from closing. Linde Noctivance® is also available as a CPAP device. Equipped with an Eco Smart™ humidifier, Linde Noctivance® produces a natural breathing environment, ensuring comfortable therapy throughout your entire sleep.

Can it be more simple?

Noctivance® is a reliable sleep device without the complicated functionality, made affordable. That's why Noctivance® keeps it simple with only one button to start therapy and one screen to check your data.

Linde is a global leader in medical gases and respiratory technology. We are confident in producing top performing, high quality, safe medical products worldwide.

Noctivance®, your solution to better health and sleep.



An affordable way to get a restful sleep

The Linde Noctivance® range offers two world-leading sleep therapy devices at affordable prices. Designed for comfort and simplicity, you can rest easy knowing these devices will help you get a better sleep.

Simple design and easy to use

- User-friendly design with large screen and buttons for ease of control and visibility
- One-click to start therapy
- Easy to assemble with simple slide-in humidifier
- Minimal effort required for cleaning

Comfortable and quiet

- Exceptionally quiet which provides rest for both you and your partner
- Ramp function increases pressure slowly, so you can comfortably fall asleep
- Adjustable humidifier prevents dry nose and throat or condensation in your mask
- Pressure relief technology that adjusts to exhalation to make breathing feel more natural



Noctivance® CPAP C100

2.4-inch display screen
Set constant pressure



Noctivance® APAP A110

3.5-inch display screen
Self adjusting pressure

5 year warranty

4–20cmH₂O pressure range to suit prescribed therapy pressure
Reslex™ exhalation pressure relief
Pressure relief technology
Mask on/off detection and power safety alerts
Clinical support with data report given to patient for progress and troubleshooting
5 year warranty for device, 6 months for consumables